



## Thoracic Catheter

Code: SR-01

Suitable for post operative drainage after Cardio-Thoracic and Thoracic surgery.

Distal end open with smooth large eyes to maximise the drainage from the chest without Trauma to the delicate lung tissues.

Proximal end fitted with specially designed tapered tongue to provide excellent forcep grip and non-traumatic insertion.

Tubes marked at 2cm from the last eye to ascertain the depth of placement.

One loose connector provided to facilitate the immediate connection of catheter with the drainage set.

Manufactured from Non-toxic Medical Grade PVC with Radio Opaque Line throughout its length.

Sterile, non-toxic & pyrogen free.

Sizes available (FG): 12,14,16,20,24,28,32,36,40

### Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box



## Thoracic Catheter with Trocar

Code: SR-02

Specially designed for quick non operative plueral and chest drainage for the relief of chest conditions such as tensions pneumothorax.

Catheter manufactured from medical grade PVC is snugly fitted over the trocar.

Blunt smooth tip.

Flat handle facilitates easy maneuverability.

Large smooth drainage eyes for efficient drainage.

Markings provided at 5, 10 ,15 & 20 cm from the distal tip help ascertain depth of penetration.

Proximal end is fitted with tapered connector for easy connection to drainage bottle.

Individually double packed in peelable pouch.