




Veterinary Infusion Sets

- Ideally suited for veterinary infusions.
- Available in vented (for bottles and bags) and non-vented (for Plastic Bottle bottles) Sharp spike for better penetration
- ISO 8536-4 compliant spike ISO 594/1 & ISO 594/2 compliant end-connection
- Options available for Latex Free
- Soft Clear Flexible Tubing
- Sterile, Non-Toxic and Pyrogen Free.
- Sterile Pack - ETO (Ethylene Oxide) Sterilization validate according to Harmonized Standards


Infusion Therapy

Available Specifications



Code: SP-01-A

101" IV Administration Set Vented
 Vented Spike 20 Drops-ml with Guard/Drip Chamber 60mm/Fluid Filter-15 micron/
 2 nos. Latex Free 'Y' Sites/Flow Regulator/Luer Tip Adaptor with Cap/3.0mm x 4.1mm
 PVC Tubing-Length 145+75+20cm/EO Sterile/Paper Poly Pack/
 Quantities per Carton = 300 pcs



Code: SP-01-B

101" IV Administration Set Non Vented
 Non-Vented Spike 10 Drops-ml with Guard/Drip Chamber 70mm/Fluid Filter-15 micron/
 2 nos. Latex Free 'Y' Sites/Flow Regulator/Luer Tip Adaptor with Cap/3.0mm x 4.1mm
 PVC Tubing-Length 145+75+20cm/EO Sterile/Paper Poly Pack/
 Quantities per Carton = 300 pcs

Veterinary Infusion Set