

# Yankaur Suction Set



Code: SR-05

Yankaur Suction Set is suitable for convenient removal of secretion, blood and debris etc per operatively.

Complete set is provided with Yankaur Suction Tip mounted on two meter long ribbed tube provided with Universal connector at both ends.

Complete set with other type of Suction tips available on request.

Universal connectors are moulded from soft PVC, so as to accommodate the Suction tip on one end and most type of connectors of Suction Apparatus at the other end.

Ribbed tube is kink resistant and has the strength to withstand the suction

## Options Available

- \* Various sizes and lengths
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box