

## Redon Drain Catheter



Code: SR-07

Soft, frosted and kink resistant PVC tubing

Markings on the tube at 19 cms from the distal tip

Provided with radio-opaque line and smooth eyes

Cross perforated length of tube : 14 cms,  
Length : 50 cms

Sizes : FG 8, 10, 12, 14, 16, 18

### Options Available

- \* Various sizes and lengths
- \* With or without XRO Line
- \* Poly Pack or Paper-Poly Pack
- \* Inner Bag or Inner Box