



# Surgical Gloves

Code: SG-01 & 02

Latex surgical gloves will help protect your hands against body fluids, help you clean up after patients and keep your work area sanitary

## Salient features:

- \* High elasticity and strength
- \* Made up of high quality natural latex
- \* Powdered or powdered free available as per requirement
- \* Moisture content: below 0.3% per glove
- \* Sterilisation by Gamma Rays / Ethylene Oxide
- \* Good design comprises bended cuffs, curved fingers for ease of use
- \* Available in sizes 6, 6.5, 7, 7.5, 8, 8

## Packing details:

Inner Packing: 1 x 50 pairs

Master Carton: 500 Pairs

Master Carton Dimension: 56 x 25.5 x 46 cm

