




EAR LOOP MASKS

CODE : MS - 17A (2-PLY)

CODE : MS - 17B (3-PLY)

- ◆ Latex Free Round Lycra Elastic 7.5" : Elongation 1:2
- ◆ Ear Loop mask is Made of Spun bond PP Fabrics
- ◆ Inside Elastic Sealing for Comfort & Neater Packing
- ◆ Easier & Quicker to wear than tie-on masks.
- ◆ International Mask Size : 175 mm X 95 mm
- ◆ Malleable PVC/PE Coated Aluminum Nose-Clip Prevents Air Leakage.
- ◆ High Bacteria / Particle / Viral Filtration Levels (> 95% , > 98%)
- ◆ Low Breathing Resistance (Delta P < 29, 4 Pa)
- ◆ Fluid Repellent Masks Available (4- ply: Splash Resistance > 120mmHg)
- ◆ Packing in PP Bags or Duplex Inner Box (Secondary Packing)
- ◆ Available in 1 / 2 / 3 / 4 ply.
- ◆ Available in Blue / Green Colors.




TIE FACE MASK

CODE : MS - 18A (2-PLY)

CODE : MS - 18B (3-PLY)

- ◆ Tie mask is Made of Spun bond PP Fabrics
- ◆ Tie mask is made with ultrasonic sealing technology.
- ◆ Tie mask is available in Blue/Green Color.
- ◆ International Mask Size: 175 mm x 95 mm
- ◆ Tie mask comes with malleable PVC/PE coated aluminum nose-clip which prevents air leakage.
- ◆ Tie bands are strong , 40 Gsm which prevents tie snapping
- ◆ Tie mask is available in 1 ply ,2 ply, 3 ply, 4 Ply.
- ◆ Packing in PP Bags or Duplex Inner Box



ACTIVATED CARBON MASK / POLLUTION MASK

Code : MS -19

- ◆ Ideal Anti Pollution Mask for commuters.
- ◆ Highly Porous 45 Gsm Active Carbon layer for excellent absorbency
- ◆ 4-Ply construction with Bacteria and Particle Filtration Efficiency > 99%
- ◆ Removes bad odors and organic vaporous
- ◆ Highly absorbent for dust and particles, yet non toxic
- ◆ Good permeability and breath ability
- ◆ Packing in PP Bags or Duplex Inner Box



VISOR MASK

Code : MS - 20

- ◆ 4-Ply Tie-On Mask with Anti Fog Optical Visor
- ◆ Black Fabric Strip over Nose-clip deflect light glare
- ◆ Long and soft PU Foam Strip under nose clip for additional user comfort and reduce fogging. Extra long nose-clip.
- ◆ Strong Tie Bands, 40 Gsm, prevents tie snapping
- ◆ 90 cm Horizontal Tie Bands accommodates all head sizes and hair lengths
- ◆ Malleable PVC coated aluminum nose-clip prevents air leakage
- ◆ High Bacteria Filtration (4-Ply: BFE>99%)
- ◆ Low Breathing Resistance (4-Ply: Delta P < 49,0 Pa)
- ◆ Fluid Repellent Masks available (4-Ply: Splash Resistance > 120 mm Hg)
- ◆ Available in Blue/Green colors.
- ◆ Application mainly for HIV Patients



**N - 95 FACE MASK /
PARTICULATE RESPIRATOR**

Code : MS - 21

- ◆ Particulates for which these respirators provide protection include contaminants such as dusts, fumes, mists and smoke.
- ◆ Filter Efficiency > 95% against 0.26 micron NaCl Aerosol airflow at 85 Liter/Minute.
- ◆ Average Maximum % Leakage of 2.5% Inhalation Resistance average of 9.4 mm H₂O.
- ◆ Exhalation Resistance average of 12.6 mm H₂O.
- ◆ 100% ultrasonically machine manufactured without manual intervention
- ◆ Excellent inhalation/exhalation performance eliminates requirement for valve and therefore potentially dangerous valve leakage
- ◆ Soft twin wired 8.5 mm aluminum nose clip molds easily on the nose to prevent leakage, mist and perspiration
- ◆ Available in different colors