

**FILTER  
ADMINISTRATION  
SETS**



- ◆ Self and fast priming of adult I.V. Set to minimize nursing time.
- ◆ Automatic air-venting in any position loss of flow due to air locks.
- ◆ Able to withstand a minimum of 45psi (3.1 bar) pressure.
- ◆ ISO 8536-4 compliant Spike & ISO 594/1 and 594/2 compliant end connection.
- ◆ Available with Flow Controller 5ml to 250ml / hr.
- ◆ Bacterial retention throughout 96 hrs. Flow.
- ◆ Options available with PVC/DEHP Free
- ◆ Filter media 1.2 / 0.2 micron
- ◆ EO Sterile, Non - Toxic and Pyrogen Free


**Available Specifications**

**Infusion Therapy**



**Code: SP-09-A**

84" IV Administration Set Vented  
Vented Spike 20 drops-ml with Guard/Drip Chamber 70mm/Fluid Filter - 15micron/Flow Regulator/  
Dial Flow Regulator/ Latex Free 'Y' Site/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/  
3.0mm x4.1mm PVC Tubing-Length 60+50+30+15cm/EO Sterile/Paper Poly Pack  
Quantities per Carton = 200 pcs



**Code: SP-09-B1**

82" IV Administration Set Non-Vented  
Non-Vented Spike 10 drops-ml with Guard/Drip Chamber 70mm/Flow Regulator/Latex Free  
'Y' Site/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/3.0mm x4.1mm  
PVC Tubing-Length 115+30+15cm/EO Sterile/Paper Poly Pack  
Quantities per Carton = 240 pcs



**Code: SP-09-B2**

82" IV Administration Set Vented  
Vented Spike 20 drops-ml with Guard/Drip Chamber 70mm/Fluid Filter - 15 micron/Flow  
Regulator/Latex Free 'Y' Site/0.2 Micron Filter/Rotating Luer Lock Male Adaptor with Cap/  
3.0mm x4.1mm PVC Tubing-Length 115+30+15cm/EO Sterile/Paper Poly Pack  
Quantities per Carton = 240 pcs



**Code: SP-09-C1**

12" Standard Bore Extension Set  
Female Luer with Cap/0.2 Micron Filter/Slide Clamp/Latex Free 'Y' Site/Rotating Luer Lock  
Male Adaptor with Cap/3.0mm x4.1mm PVC Tubing-Length 10+10+10cm/EO Sterile/  
Paper Poly Pack / Quantities per Carton = 300 pcs



**Code: SP-09-C2**

12" Standard Bore Extension Set  
Female Luer with Cap/1.2 Micron Filter/Slide Clamp/Latex Free 'Y' Site/Rotating Luer Lock  
Male Adaptor with Cap/3.0mm x4.1mm PVC Tubing-Length 10+10+10cm/EO Sterile/  
Paper Poly Pack / Quantities per Carton = 300 pcs

**Extension Lines**