

Nasopharyngeal Airway



Code : AS-21

- ◆ Nasopharyngeal Airway can be used on breathing or non-breathing patients. However care must be taken to establish the level of consciousness as a gag reflex may or may not be present.
- ◆ Suitable for suctioning in I.C.U.s & provide an alternative way when oropharyngeal airway is occluded.
- ◆ Made of Medical Grade PVC readily conforms to the nasopharyngeal anatomy.
- ◆ Thin wall construction provides maximum inner diameter.
- ◆ Interior walls with channel grooves, facilitate easy insertion of catheter.
- ◆ Specially designed bevel tip reduces mucosal trauma during insertion.
- ◆ Available Sizes : 2.5 , 3.0 , 3.5 , 4.0 , 4.5 , 5.0 , 5.5 , 6.0 , 6.5 , 7.0 , 7.5 , 8.0 , 8.5
- ◆ Individually packed in Paper Poly Pack.